

## A New Gold Standard in Infusion Management

# niki T34 Ambulatory Syringe Pump



The T34 Syringe driver is the new standard in ambulatory Syringe drivers. Offering a comprehensive range of features designed to maximize safe, economic and simple administration of medication via multiple program options.

The T34 is light, compact & robust ensuring it can be used in a wide variety of treatment settings including hospital in- and out-patient & homecare.

The T34 constitutes a major advance in technology incorporating active & passive safety features configurable to local practice and protocols and lockable to prevent use or change by unauthorized personnel.

These features include security access codes, free-flow protection (optional), downloadable patient and event history and automatic detection of a wide range of syringes from 2 to 50 ml. The T34 software can be adjusted to local users needs including choice of language for information display.

Manufactured by: Caesarea Medical Electronics Limited  
Address: Staufenburgstr. 23, Lichtenstein  
POB 1248, Lichtenstein 72805, Germany.

©2004 by Caesarea Medical Electronics Ltd. All rights reserved.

[www.cme-infusion.com](http://www.cme-infusion.com)

FDA Approved 510(K) : K050800

- **Lightweight & Compact** – The T34 is the most compact & lightest pump in its class.
- **Optional pole-mounting in a docking station**
- **Comprehensive range of pump accessories** including shoulder carry bag, lockable syringe cover and docking station.
- **3 levels of Access Security Codes**, ensures only authorised personnel can use the device and only senior personnel can set up the devices parameters.
- **Downloadable 24 Hour Patient History** (including charts) & **1000 Event Log** for full accountability & technical data.

CE 0473

# Niki T34 Ambulatory Syringe Driver

Rate:	From 0.1 to 650 ml/h
Time Range	From 10 Minutes to 100 Hours
Dimensions	165 x 23 x 40mm (6.5 x 0.9 x 1.8 inches)
Weight	247 grams including batteries
Battery Type	Alkaline 9v Batteries
Battery Life	100 full Syringes
Controls	Graphic Display 122 x 32 pixels with back light (8 keys) & dual color operation LED
Syringe range	2ml to 50 ml syringe (any brand)
Languages	Selectable- English, Italian, Spanish

## Alarms

Occlusion	Downstream (between pump & patient)
End Program	Program Completed
Near End	Program near completion
End Travel	Near End Syringe, Empty Syringe
Low Battery	Battery charge is lower than 7v
End Battery	Battery Depleted
Check program	Wrong Setting - (protocol requested is outside of safe limits set by senior personnel under level three access code)
High Pressure	The pressure in the line has exceeded the pre-set limit

## Features

Automatic detection of the syringe size and brand

Programs : Continuous (Dose – Optional)

Serial Port (RS232) connection (Optional)

