



WARMFLO* 588

Blood and Fluid Warmer



*Fluid Warming
has never been easier*

tyco

Healthcare

Warmflo* 588

Blood and Fluid Warmer



FLUID WARMING HAS NEVER BEEN MORE PREDICTABLE

- Ideal for all routine clinical situations at a wide range of temperature settings; up to 42°C
- New WF-100 fluid warming cassette design provides clean and reliable warming

FLUID WARMING HAS NEVER BEEN EASIER

- Simply slide in the disposable fluid warming cassette and the system is operational
- Compact and portable - can be mounted on existing IV poles
- Versatile - WF-100 fluid warming cassettes work in Warmflo* 588 and 537 fluid warmers

AUTOMATIC PRESSURIZED INFUSION CHAMBER

This fully automatic chamber allows for simple, reliable and fast infusion. The pressure level is maintained with a unidirectional valve whenever the chamber is being moved. Cleaning is facilitated and bacteria propagation is avoided with smooth streamlined surfaces. The chamber is available in 500 ml and 1000 ml versions; it is equipped with a pressure gauge / governor as well as an ON/OFF switch.

FEATURES AND BENEFITS:

- Provides reliable, low-maintenance, dry-heat technology
- Advanced system design offers easy-to-use fluid warming cassette
- Fully compatible with in-line IV hand pumps

- *Infusion chamber operating under AD Merkblätter pressure (German standard) with compressed air source (wall mounted or cylinder)*

REF	Description	Quantity
AD 500	Quick connection / disconnection	1
AD 1000	Quick connection / disconnection	1

- *Fluid and disposable items heater*

REF	Description	Quantity
FW588	Fluid Warmer	1
WF-100	Fluid Warming Cassette	20
WF-250	Fluid heating cassette equipped with blood filter	20

*Trademark © 2003 Tyco Healthcare Group LP. All rights reserved.

C-TM-WF 588/GB



Tyco Healthcare UK Limited
154 Fareham Road
GOSPORT Hampshire
PO 13 OAS
UK
Tel: +44 (0) 1329 224226
Fax: +44 (0) 1329 224334

Tyco Healthcare
Regional Export
154 Fareham Road
GOSPORT Hampshire
PO 13 OAS
UK
Tel: +44 (0) 1329 224258
Fax: +44 (0) 1329 224086

Tyco Healthcare Nederland BV
Hogeweg 105
NL-5301 LL ZALTBOMMEL
The NETHERLANDS
Tel: +31 418 57 66 00
Fax: +31 418 57 67 91
TycoHealthNed@TycoHealth.com
www.tycohealthcare.nl

Tyco Healthcare
Birchwood Court Montrose Road
P.O. Box 8108, Halfway House
1685 MIDRAND
SOUTH AFRICA
Tel: +27 (0) 11 265 2000
Fax: +27 (0) 11 265 2050

Tyco Healthcare Belgium N.V.-S.A.
Generaal de Wittelaan 9/5
B-2800 MECHELEN
BELGIUM
Tel: +32 (0) 15 29 44 50
Fax: +32 (0) 15 29 44 55
tycohealthcare.Belgium@tyco-healthcare.com

Tyco Healthcare Norden AB
Gårdsvägen 18, Box 710
SE-169 27 Solna
SWEDEN
Tel: +46-8 585 605 00
Kundservice +46-8 585 605 05
Fax: +46-8 585 605 06

Tyco Healthcare Danmark
Langebrogade 6E, 4.th.
DK-1411 København K
DENMARK
Tel: +45-70 27 52 50
Kundeservice +45-70 27 52 50
Fax: +45-70 27 51 50

Tyco Healthcare Norge
Løkketangen 20
N-1337 Sandvika
NORWAY
Tel: +47-67 80 55 10
Kundeservice +47-67 80 55 18
Fax: +47-67 80 55 11

Tyco Healthcare Finland
Köydenpunojankatu 2D, 6 krs, PL340
FIN-00181 Helsinki
FINLAND
Puh +358-(0)9-6226 840
Asiakaspalvelu
+358-(0)9-6226 8418
Fax: +358-(0)9-6226 8460